Geophysical Research Abstracts, Vol. 9, 05702, 2007 SRef-ID: 1607-7962/gra/EGU2007-A-05702 © European Geosciences Union 2007



Fluid inclusion hydrogen and oxygen isotope analyses using the "Amsterdam Device": a progress report

H.B. Vonhof (1), T.C. Atkinson (2), M.R. van Breukelen (1) and O. Postma (1)

(1) Vrije Universiteit Amsterdam, Faculty of Earth and Life Sciences, Amsterdam, Netherlands (2) University College London, Bloomsbury Environmental Isotope Facility, London, U.K. (hubert.vonhof@falw.vu.nl)

The Amsterdam Device is a continuous-flow preparation device for on-line δ^2 H and δ^{18} O analysis of fluid inclusions in speleothem carbonate. The design is based on a relatively low-cost adaptation to a Thermo Finnigan TC-EA pyrolysis furnace. Standard specifications of the Thermo Finnigan TC-EA allow for simultaneous δ^2 H and δ^{18} O analysis of 0.2 microliter of water, reproducible within 2 permil for δ^2 H and 0.5 permil for δ^{18} O (1SD). The technical description of the Amsterdam device was published in 2006, together with the results of some first crushing experiments. These experiments showed that the technique works well for δ^2 H analyses, but not enough experiments were done at the time to evaluate the applicability of the device for δ^{18} O analyses of the same samples.

Over the last year we have run extensive experiments with the so-called jump routine, which allows combined $\delta^2 H$ and $\delta^{18} O$ analysis of the water released from a single crush. The jump routine is a standard feature in the software for our Thermofinnigan TC-EA and Delta XP, and involves a calibrated rapid magnet jump of the mass spectrometer in the time between the entry of the GC-separated H₂ and CO peaks.

In these recent experiments we have made two main changes compared to the Vonhof et al (2006) protocol:

1: We have installed a reverse-flow carrier gas setup, which has decreased the memory effect on TC-EA analyses (Gehre et al., 2004). To accommodate that, we split the incoming carrier gas in two flows: A) a ~ 10 ml/min flow entering at the bottom connector of the TCEA reactor, and B) a ~ 100 ml/min carrier flow that leads through the crusher and cold trap to enter in the top of the glassy carbon reactor.

2: For the freezing protocol we apply an ethanol slush at ~ -90 degrees Celsius, instead of liquid nitrogen applied in the older experiments. This is essential because isotopically distinct CO₂ released from the crushed CaCO₃ sample is trapped along with the inclusion water under liquid nitrogen temperature and contributes to the δ^{18} O value analysed with the pyrolysis method. At ethanol slush temperatures, water will be trapped quantitatively and CO₂ will not, leading to a inclusion water δ^{18} O analysis unaffected by CO₂ derived oxygen.

With this new protocol we have analysed speleothem material from different settings and a selection of results is presented. Often, $\delta^2 H$ and $\delta^{18} O$ data plot near the Global Meteoric Water Line (GMWL), which lends support to the accuracy of the analyses. If data plot away from the GMWL it is usually the $\delta^{18} O$ data that are off-set. These $\delta^{18} O$ off-sets are believed to be analytical artefacts. Our experiments did not provide evidence for oxygen exchange between fluid-inclusion water and surrounding carbonate.

References:

Gehre, M., H. Geilmann, J. Richter, R. A. Werner, and W. A. Brand, 2004, Continuous flow ${}^{2}\text{H}/{}^{1}\text{H}$ and ${}^{18}\text{O}/{}^{16}\text{O}$ analysis of water samples with dual inlet precision: Rapid Communications in Mass Spectrometry, v.18, p. 2650–2660

Vonhof, Hubert.B., Martin van Breukelen, Onno Postma, Peter J. Rowe, Tim C. Atkinson and Dick Kroon, 2006, A continuous-flow crushing device for on-line δ^2 H analysis of fluid inclusion water in speleothems: Rapid Communications in Mass Spectrometry, v. 20, p. 2553-2558.